


Radiology Image Share™ Smart Solutions

Virtua® Smart Media™



Share your DICOM studies
via the cloud, simply & securely



CODONICS

Defining Image Distribution

Radiology Image Share™ Smart Solutions

Virtua® Smart Media™

Share DICOM Studies via the Cloud, Simply and Securely

Virtua Smart Media provides your facility with a simple way to present your DICOM studies through a virtual medium (Cloud technology). This streamlined and secure service enhances workflow efficiency and optimizes data sharing throughout your hospital.

With Virtua Smart Media, the Cloud is Yours

You select the studies/series you would like to send and then enter the recipient's e-mail address (manually or from a directory). An e-mail is sent automatically, providing the recipient with access to all images in a simple click. Images are available via the user-friendly viewer (compatible with all PCs, without software installation).



Key Features

- Sends e-mails to multiple recipients
- Enables easy contact management
- Keeps patient information anonymous
- E-mail template (HTML) is fully configurable
- User selects the specific series for the study
- Fully configurable downloads are provided with an expiration period

Completely Secure

Thanks to the SS protocol and encrypted data (Advanced Encryption Standard 256), Virtua Smart Media provides optimum security. The storage space is certified SSAE16/SOC1, SOC2, ISAE 3402. In the standard version, Virtua Smart Media enables you to store 10 studies simultaneously in the Cloud – when you send the eleventh study, the first one is deleted in accordance with the FIFO (First In, First Out) methodology.

Hardware sourced separately. Please speak to your Codonics sales representative for more details.



Call Codonics today at **+1.440.243.1198** or visit us at www.codonics.ae, www.codonics.lat, or www.codonics.fr for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. www.codonics.com/ip/patents/
Copyright © 2024-2025. Codonics Inc. 1/2025



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com

900-904-002.01