


Radiology Image Share™ Smart Solutions

DP Smart™ (DICOM Print)



CODONICS

Defining Image Distribution

Radiology Image Share™ Smart Solutions

DP Smart™ (DICOM Print)


DP Smart

DP Smart enables you to print studies from any DICOM modality to any paper printer powered by Windows. The solution makes it easy to output images directly to paper, replacing expensive film. Now, with DP Smart, your institution can provide hard copy study reference to your patients and referring physicians quickly and cost-effectively.

Key Features

- Compatible with any printer on the market
- Convert DICOM printing to Windows printing
- Print studies from any DICOM modality
- Multi-format printing with multiple DICOM Print AE Titles: A4, A3, color, black and white, etc.
- Customizable header and footer (adding/editing text, logo, etc.)
- Print adjustment panel according to the connected printer

Hardware sourced separately. Please speak to your Codonics sales representative for more details.



Call Codonics today at **+1.440.243.1198** or visit us at www.codonics.ae, www.codonics.lat, or www.codonics.fr for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. www.codonics.com/ip/patents/
Copyright © 2024-2025. Codonics Inc. 1/2025



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com