

Radiology Image Share™ Smart Solutions

P2P Smart™ (Print to PACS)



Import documents
into PACS in Just ONE Click



CODONICS

Defining Image Distribution

Radiology Image Share™ Smart Solutions

P2P Smart™ (Print to PACS)

Import Documents into PACS in Just ONE Click

Print to PACS is the ideal solution for importing and integrating documents into your PACS archiving system. Easy to use, the Print to PACS solution converts your printed documents into DICOM images and assigns them to a new study, a stored study on a DICOM server, or a patient study retrieved from a DICOM Worklist.

Key Features



- Convert Windows printing to DICOM images
- Quickly send Windows printable documents to PACS
- Select, edit and sort images before sending to PACS
- Create a study based on a DICOM Worklist entry or an existing patient in the PACS
- Assign to an existing study in the PACS
- Configure easily and simply
- Manage users
- View sent history

Connectivity and DICOM Services

Fully DICOM 3.0 compliant and provides the following DICOM services:

- C-ECHO SCU
- C-FIND SCU
- C-FIND MWL SCU
- C-STORE SCU

Hardware sourced separately. Please speak to your Codonics sales representative for more details.



Call Codonics today at **+1.440.243.1198** or visit us at www.codonics.ae, www.codonics.lat, or www.codonics.fr for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. www.codonics.com/ip/patents/ Copyright © 2024-2025. Codonics Inc. 1/2025



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
Email info@codonics.com
www.codonics.com